USER MANUAL

CONTENTS

Guarantee	2
Introduction	4
Identification data	6
Component List	6
Installation	7
Controls	8
Cleaning, Disinfection & Sterilization	10
Regular maintenance	13
Replacement of the fuses	13
Delivery	14
Disposal procedures and precautions	14
Tips	14
Symbol description	15
Trouble shooting	15
Technical data	15
After service	15
Environmental protection	15
Manufacturer's right	15

GUARANTEE

- 1) Before being placed on the market, all equipment undergoes a thorough final check to ensure that it is in proper working order.
- 2) We guarantee its products, purchased new from a dealer or importer, to be free from manufacturing or material

defects for:

- TWO YEAR from the date of purchase for the device;**
- ONE YEAR from the date of purchase for the handpiece with its cord.**

- 3) Throughout the warranty period, undertakes to repair (or, at their sole discretion, to replace) free of charge any parts that, in their opinion, are faulty.

Complete replacement of products is excluded

- 4) Cannot accept any liability for direct or incidental damage or personal injury in the following cases:
 - a) If the equipment is used for purposes other than that for which it is intended;
 - b) If the equipment is not used in accordance with all the instruction and requirement described in this manual;
 - c) If the wiring system in the room where the equipment is used does not comply with the applicable standards and appropriate requirements;
 - d) If any assemble operations, extensions, settings, alterations or repairs have been carried out by personnel not authorized;

- e) If the environmental conditions in which the device is kept and stored do not comply with the requirements indicated in the chapter on technical specifications.

- 5) Accidental damages due to transport, incorrect use or carelessness or to connection to power supplies other than as envisaged and damage to the signaling lamps handpiece and all accessories are excluded from the warranty. The warranty will no longer apply if the apparatus has been tampered with or repaired by unauthorized personnel.

- 6) **Warning:**

The warranty is valid only if the warranty slip enclosed with the product has been completed in full and returned to us or, if appropriate, to your dealer or importer within 20 days from the date of purchase, as proven by the consignment note/invoice issued by the dealer/importer.

In order to benefit from the warranty service, the customer must return the apparatus to be repaired to the dealer/importer from which it was purchased, at his own expense.

- 7) The apparatus should be returned suitable packed (possibly in its original packing material).

- 8) Accompanied by all the accessories and by the following information:

- a) Owner's details, including his telephone number;
- b) Details of the dealer/importer;
- c) Photocopy of the consignment note/purchase invoice of the apparatus issued to the owner and indicating, in addition to the date, also the name of the apparatus and its serial number;
- d) A description of the problem.

- 9) Transport and any damages caused during transport are not covered by the warranty. In the event of failure due to accidents or improper use, or if the warranty has lapsed, repairs to produce will be charged on the basis of the actual cost of the materials and labour required for such repairs.

AI SURGERY PRO WARRANTY CARD			AI SURGERY PRO WARRANTY CARD			AI SURGERY PRO WARRANTY CARD		
Name of Customer			Name of Customer			Name of Customer		
Address Details			Address Details			Address Details		
Postal Code			Postal Code			Postal Code		
Tel			Tel			Tel		
Model			Model			Model		
Unit No.			Unit No.			Unit No.		
Handpiece No.			Handpiece No.			Handpiece No.		
Purchase Date			Purchase Date			Purchase Date		
Contact person			Contact person			Contact person		
Date	Maintenance Record	Repairer	Date	Maintenance Record	Repairer	Date	Maintenance Record	Repairer
1.For Customer 2.For Distributor 3.Return to Manufacturer			1.For Customer 2.For Distributor 3.Return to Manufacturer			1.For Customer 2.For Distributor 3.Return to Manufacturer		

INTRODUCTION

1.1. Foreword

Before proceeding with the installation, use, maintenance or any other activities on the equipment please read the manual carefully.
Important: To avoid causing personal injuries or damages to property, read all the points concerning "safety requirement" contained in this manual with particular attention.

Depending on the level of risk involved, safety requirements are classed under the following indications:

Danger(always referred to personal injury)

Warning(referred to possible damage to property)

Not to position the device to make it difficult to operate the disconnection device. In the presence of electromagnetic interference environment, the device may be malfunctioning. Do not install AI Surgery Pro near equipment that releases magnetic waves.

Device with electromagnetic launcher will affect the normal operation of AI Surgery Pro, do not run both devices at the same time.

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

The purpose of this manual is to ensure that operators are aware of the safety requirements, of the installation procedures and of the instructions for correct use and maintenance of the apparatus.

The user is not authorized to tamper with the equipment under any circumstances.

If any problems are encountered, please contact a Manufacturer Service Centre.

Any attempts on the part of the user or any unauthorized personnel to tamper with or alter the apparatus will invalidate the warranty and release the Manufacturers from any liability in respect of any harm or damage to persons or property.

The information and illustration contained in this manual are up-dated to the date of publication indicated on the last page.

We are committed to continuous up-dating of the products, which may entail changes to components of the equipment. If there are any discrepancies between the descriptions contained in this manual and your equipment, please contact your dealer or the Manufacturer after-sale service for explanations.

Using this manual for purposes other than those relating to the installation, use and maintenance of the equipment is strictly prohibited.

1.2. Description of the Device

Thanks to its controlled three-dimensional ultrasound oscillations, the original AI Surgery Pro technique rings in a new age for osteotomy and osteoplasty in Implantology, Periodontology, Endodontics, and Orthodontic Surgery.

It's main features are:

Micrometric cutting: Maximum surgical precision and intra-operative sensibility;

Selective cutting: Maximum safety for the soft tissues;

Cavitation effect: Maximum intra-operative visibility (bloodless field);

The equipment has an automatic tuning circuit that offsets wear of the tips, thus ensuring work in constant conditions of maximum efficiency.

1.3. Intended Use

The AI Surgery Pro is a piezoelectric device for bone surgery that enables osteotomy and osteoplasty techniques to be applied to in almost any anatomical situation. This equipment can be used in the following fields:

- | | | |
|----------------------|------------------------|---------------------------|
| a) Oral surgery; | b) Orthopedic surgery; | c) Maxillofacial surgery; |
| d) Cosmetic surgery; | e) Neurosurgery; | f) Otolaryngology. |

This equipment cannot function in places where there is an inflammable atmosphere (anaesthetic mixture, oxygen, etc).

1.4. Safety requirements

We will not accept any liability for direct or incidental personal injury or damage to property in the following cases:

- 1) If the equipment is used for purposes other than that for which it is intended;
- 2) If the equipment is not used in accordance with all the instructions and requirements described in this manual;
- 3) If the wiring system in the room where the equipment is used does not comply with the applicable standard and appropriate requirements;
- 4) If any assembly operations, extensions, settings, alterations, or repairs have been carried out by personnel not authorized;
- 5) If the environmental conditions in which the device is kept and stored do not comply with the requirements indicated in the chapter on technical specifications.

Danger: Qualified and specialized personnel.

This equipment may be used only by specialized and suitably trained personnel such as surgeons. If correctly used, this equipment does not give rise to side effects. Improper use, on the other hand, will give rise to the transmission of heat to the tissues.

Danger: Intended use.

Use the equipment solely for the purpose for which it is intended (see point 1.3), failure to comply with this requirement could lead to serious harm to the patient and/or to the operator and/or damage to/failure of the equipment.

Danger: Contraindications.

Do not use the AI Surgery Pro on patients with pace-makers or other implantable electronic devices. The same requirement applies also to the operator.

Danger: Contraindications.

An electrosurgical knife could interfere with the correct functioning of the device.

Danger: Cleaning, disinfection and sterilization of new or repaired products.

All new or repaired products are delivered in no sterile conditions. Before being used for treatments, all new or repaired products should be cleaned, disinfected and sterilization following the instructions provided under point 8 strictly.

Danger: Use only original accessories and spare parts.

Danger: Check the condition of the device before treatment.

Always make sure that there is no water under the apparatus. Before each treatment always check that the equipment is in proper working order and that the accessories are efficient. Do not carry out the treatment if any problems are encountered in operating the device. If the problems concerning the equipment contact an authorized technical service center.

Danger: Breakage and wear of the tips.

The high-frequency vibrations and wear may, very occasionally, lead to breakage of the tip. Tips of which the shape has been changed or which are otherwise damaged are liable to break during use. Any such tips should definitely not be used. It is necessary to instruct the patient to breathe through his nose during the treatment in order to avoid ingestion of the broken off fragment of the tip.

Danger: Do not install this equipment anywhere there is a risk of explosions.

This equipment cannot function in places where there is an inflammable atmosphere. (anaesthetic mixture, oxygen, etc)

Danger: Personnel injury.

The foot switch of the AI Surgery Pro must not be activated when the door of the peristaltic pump open. (Fig.3-Ref.A).Moving parts could injure the operator.

Danger: Contraindication.

Do not carry out this treatment on metal or ceramic prosthetic artifacts. The ultrasonic vibrations could lead to decrementing of such artifacts.

Danger: Contraindication.

After autoclave sterilizing of the handpiece, wait for it to cool down completely before using it.

IDENTIFICATION DATA

COMPONENTS LIST

2.1. Identification data

An exact description of the model including the serial number of the equipment will make it easier for our After-Sale Service to respond quickly and efficiently to your enquiry.
Always provide the above information whenever you contact a Manufacturer Service Centre.

2.2. Data plate of the device

Each device has its own data plate (Fig.1), on which technical specifications and serial numbers are indicated. The data plate is on the rear of the device.

The remaining data are included in this manual (see Page 16).


Fig.1

2.3. Data plate of the handpiece

The serial number of the AI Surgery Pro handpiece is engraved on the ring nut (Fig.2)

Fig.2

3. The materials included in the supply may vary in the case of promotional campaigns.

Device	Wireless Foot Switch	Foot Switch Power Adapter	Tubes
Water bottle Bracket	Handpiece with Cable	Handpiece holder	Connector
Torque wrench	Sterilize box	Tip Holder and Tips	Power-supply cable

INSTALLATION

4.1. Safety requirements during Installation

Danger: The wiring system of the premises where the apparatus is installed and used must comply with the applicable standards and the relevant electrical safety requirements.

Danger: Do not install the apparatus in places where there is a risk of explosion. The apparatus may not be used in areas where there are inflammable atmospheres (anaesthetic mixtures, oxygen, etc).

Danger: Install the apparatus in a place where it will be protected from blows and from accidental sprays of water or other liquids.

Danger: Do not install the device on or in the vicinity of sources of heat. Install it such a way that there is an adequate circulation of air around it. Leave sufficient free space around it, in particular with reference to the fan on the rear. (Fig.4)

Warning: Do not expose the apparatus to direct sunlight or to sources of UV light.

Warning: The apparatus is transportable, however it must be handled with care when it is moved.

Warning: Before connecting the cord to the device, make sure that the electrical contacts are perfectly dry. If necessary, dry them with the air syringe.

Warning: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Warning: The capacity of brine bottle hang on brine bottle holder should not be more than 1000ml, and the weight should not exceed 1Kg.

4.2. Initial installation

To ensure perfect operation of the equipment, it is installed by technical personnel authorized. The equipment will be installed in a suitable and handy place for it to be used. The technician must:

- 1) Install the device in a suitable place;
- 2) Explain the main aspects of correct installation to the user;
- 3) Fill in the installation form, including the purchaser's data;

4.3. Connection the accessories

The accessories listed as follow should be connected with the AI Surgery Pro:

- 1) Insert the silicone tube into the peristaltic pump, proceeding as follows:
 - a) Open the door(Fig.3)as far as it will go.
 - b) Position the tube in the impeller.
 - c) Close the door completely.

Danger: personnel injury.
The footswitch of the AI Surgery Pro must not be activated when the door of the peristaltic pump open .Moving parts could injure theoperator.

- 2) Insert the rod for supporting the bag into the holes provided for it (Fig.4-Ref.A);
- 3) Plug the power cable into the connector on the casting of the device (Fig.4-Ref.D) and then into the power outlet;
- 4) Insert the tube of AI Surgery Pro cord to the cord connector on the device (Fig.4-Ref.B);
- 5) Connect end of the tube of the peristaltic pump;
- 6) Connect the flow-control system to the bag containing the appropriate liquid for the treatment;

7) Use the torque wrench to screw the tip (Fig.5) till the clattering voice;

Fig.4

8) Press the button "on/off" (Fig.4—Ref.C), then can use the device.


Fig.5

4.4. Testing of the device

All the devices are checked and tested by , including all the parts. When testing, all the parts will work in intermittent operation. The test emphasized that all the problems are from the failure parts. This procedure ensures the function and reliability of all the parts.

CONTROLS

5.1. Description of the controls

This section illustrates the parts of the front panel of the AI Surgery Pro unit, enabling the controls described in this manual to be located immediately.


Fig.6

1	Bone function	2	Endo function
3	Perio function	4	Custom function
5	Clean function	6	Display the status of foot switch
7	Display handpiece connection status	8	Display Water flow rate, press to adjust the water flow rate.
9	Decrease, adjust the Power	10	Increase, adjust the Power
11	Display the current Power rate		

In Clean function, press the foot switch, the device can clean the tube *(Recommended at least 30 seconds)*

5.2. Safety requirements during use

Danger: Contraindications.

Do not use the AI Surgery Pro on patients with pacemakers or other implantable electronic devices. This requirement also applies to the operator.

Danger: Breakage and wear of the tips.

The high-frequency vibrations and wear may, very occasionally, lead to breakage of the tip. Tips of which the shape has been changed or which are otherwise damaged are liable to break during use. Any such tips should definitely not be used. It is necessary to instruct the patient to breathe through his nose during the treatment in order to avoid ingestion of the broken fragment of the tip.

Danger: Control of infections.

For maximum safety of both the patient and the operator, clean, disinfect and sterilize the piezo electronic handpiece, the tips and the torque wrench after each treatment.

Warning: Contraindication.

Do not carry out this treatment on metal or ceramic prosthetic artifacts. The ultrasonic vibrations could lead to decremementing of such artifacts.

Warning: Contraindication.

After autoclave sterilizing of the handpiece, wait for it to cool down completely before using it.

Warning: The electrical contacts inside the cord connector must be dry.

Before connecting the handpiece to the device, make sure the electrical contacts of the connector are perfectly dry, in particular after the autoclave sterilization cycle. If necessary, dry the contacts by blowing air onto them with the syringe.

Warning:

To use the device correctly, it is necessary to press the footswitch and start it up without letting the tip rest on the part to be treated. This will allow the electronic circuit to detect the point where resonance of the tip is without any interference, thus enabling optimum performance.

If this is not done, contact with the part to be treated or with other surfaces before start-up could cause tripping of the protection systems.

Warning: For spray treatment, use only tips through which liquid is passed.

5.3. Instruction for use

- 1) Open the air intake on the drip system;
 - 2) Screw the chosen tip onto the AI Surgery Pro handpiece until it is flush against it;
 - 3) To use the torque wrench correctly (Fig.5) proceed as follow;
 - a) Hold the body of the handpiece firmly;**Warning:**
Do not grip the end part of the handpiece or the cord, only the plastic casting (Fig.5) and do not turn it while fastening the tip in place;
 - b) Turn the wrench in a clockwise direction until the cultch engages (till making clicking sound);
 - c) The tip is now properly tightened in place;
- 4) Make sure that the AI Surgery Pro handpiece is correctly connected to the handpiece connector (Fig.4-Ref.B);
 - 5) Check the display to see the power level that has been set, if the type of power required differs from the level that has been set, use the key "+" and "-" on Multi-function foot pedal for selecting, depending on the type of function that has been set;
 - 6) Check the display to see the delivery rate of the peristaltic pump, if the delivery rate required is other than the level that has been set, use the key "Water" on Multi-function foot pedal to choose, depending on the type of function that has been set.

5.4. Rules for keeping the device in proper working order

- 1) Check the state of wear of the tips periodically and replace any for which a drop in performance is noted;
- 2) Do not alter the shape of the tips by bending or filling them;
- 3) Replace any tip that has become deformed or damaged by impacts;
- 4) Always make sure that any threaded parts and their contact surfaces are perfectly clean;
- 5) If a tip becomes too worn, the device will stop working.

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